



Arm Blood Pressure Monitor

Model:U807



Instruction Manual

Table of Contents

Introduction	3
Safety Information	4
Product Structure	
--Each part name	7
Battery installation	8
Setting mode	
-- How to set	10
Proper use of the unit	
--Pre-measurement	13
--Common factors	
of wrong measurement	13
--Fitting the cuff	14
--Measuring procedure	15
--Discontinuing a measurement	15
--Memory-recall of measurement	15
--Read memory record	15
--Memory-clear of measurements	16
--Checking IMEI details	16
--Signal Strength Indicator	17
About blood pressure	17
Exceptional situations	19
Care and maintenance	20
Specification	21
Warranty information	22
EMC Declaration	23
FCC Statement	27

2

Introduction

▲ Your new digital blood pressure monitor uses the oscillometric method of blood pressure measurement. This means the monitor detects your blood's movement through your brachial artery and converts the movements into a digital reading. An oscillometric monitor does not need a stethoscope, so the monitor is simple to use.

▲ Intelligent inflation will reduce the uncomfortable feeling by incorrect inflation, and shorten the measurement time, prolong the cuff's usage lifetime.

▲ 2x90 sets memory function, each measurement result will be displayed on the screen, and automatically stored. This unit has blood classification index, could easy to check your classification index, could easy to check your classification index.

▲ Please read the manual carefully before you use the unit, and keep the manual well after using.

CONTRAINDICATION

This product can't be used in patients who is with severe heart insufficiency to avoid suffocation and death. This product is not suitable for infants and children.

INTENDED USE

This automatic blood pressure monitor intends to measure the systolic pressure, diastolic pressure and pulse rate through upper arm. It's expected to be used at home or in the hospital, intended for people over 12 years old.

3

Safety Information

■ To assure the correct use of the product, basic safety measures should always be followed including the warning and the caution listed in the instruction manual:

Symbol descriptions

The following symbols may appear in this manual, on the label, on the device, or on its accessories. Some of the symbols represent standards and compliances associated with the device and its use.

⚠ **WARNING:** This alert identifies hazards that may cause serious personal injury or death.

⚠ **CAUTION:** This alert identifies hazards that may cause minor personal injury, product damage or property damage.

⚠ **Type BF applied part**

🏭 **Manufacturer**

SN **Specifies serial number**

⚡ **Authorized Representative in the European Community**

♻ **DISPOSAL:** Do not dispose of this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.

⚡ **Direct current**

📖 **Follow instructions for use**

4

Safety Information

⚠ Those who have arrhythmia, diabetes, blood circulation or apoplexy problem, please use under the physician's instruction.

⚠ Contact your physician for specific information about your blood pressure. Self diagnosis and treatment which use measured results may be dangerous. Follow the instructions of your physician or licensed healthcare provider.

⚠ Please place on a high place where children can't be touched.

⚠ No modification of this equipment is allowed.

⚠ Do not modify this equipment without authorization of the manufacturer.

⚠ If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of equipment.

⚠ The cuff hose around neck may cause the suffocation.

⚠ The swallowing of small part like packaging bag, battery, battery cover and so on may cause the suffocation.

⚠ Please don't use a dilution agent, alcohol or petrol to clean the unit. Please don't hit heavily or fall down the product from a high place. Use the right cuff, otherwise it can not work.

⚠ Do not replace or remove the battery from device (in the case of device with rechargeable lithium battery).

⚠ Do not use a cellular phone near the unit. It may result in operational failure.

⚠ Please avoid using in high radiant area in order to make your measuring data correctly.

5

Safety Information

⚠ Do not use the equipment where flammable gas (such as anesthetic gas, oxygen or hydrogen) or flammable liquid (such as alcohol) are present.

⚠ **WARNING:**

♻ Do not dispose of electrical appliances as unsorted municipal waste. Use separate collection facilities. Contact your local government for information regarding the collection systems available. If electrical appliances are disposed of in landfills or dumps, hazardous substances can leak into the groundwater and get into the food chain, damaging your health and well-being.


Classification

- Internally powered equipment;
- Type BF applied part;
- Protection against ingress of water or particulate matter: IP21;
- Not category AP / APG equipment;
- Mode of operation: Intermittent operation.

⚠ The user must check that the equipment functions safely and see that it is in proper working condition before being used.

6

Product Structure



7

Battery installation

SIM card installation

Remove the battery cover from the battery compartment, the SIM card slot is on the right side (as picture shown). Insert SIM card into the slot in accordance with direction " " and chip side down.

Battery installation

The device is equipped with 1pc 3.7V 1000mAh rechargeable lithium battery.

Noted:

⚡ Do not replace or remove the battery from device.

LOW battery and charge

When power on, the low battery symbol will be displayed once the device starts. Plug in the device using the included USB cable to charge.

Take the USB cable and connect it to the charger (charge connector) or to your USB port on your computer or power bank. Connect the end of the round tip to the USB socket of device.

Noted: the lithium battery could be used

for about 70 times when charge is full. Charge time is 4-6 hours.

During charging when power off, the LCD will shows the blood pressure bars on the left side from one bar to six bars in a continuous loop. And unmoved 6 bars along with the battery symbol will be displayed once battery is full, as picture shown.



8

Battery installation

If the device under setting mode and battery charging in the same time, LCD can show the same as above, e.g. user setting



⚠ **WARNING:**

Dispose of the battery in accordance with all federal, state and local laws. To avoid fire and explosion hazard, do not burn or incinerate the battery.

Adapter Usage

- When optional AC adapter is in use, it should comply with the requirement of IEC 60601-1:2005. Furthermore, all configurations shall comply with the requirements for medical electrical systems (see IEC 60601-1-1 or clause 16 of the 3Ed. of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical equipment configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems. Attention is drawn to the fact that local laws take priority over the above mentioned requirements. If in doubt, consult your local representative or the technical service department.
- When using AC power, to avoid possible damage to the monitor, use only the exclusive AC adapter that can be purchased from authorized dealers. Other adapters may vary in output voltage and polarities.
- Insert the adapter plug into the hole on the backside of the unit as pictured

9

Setting mode

- Insert the other side of the adapter into the outlet with 100-240V.
- To remove the AC adapter, disconnect the adapter plug from the outlet first and then disconnect the cord from the unit's socket.

Adapter technical features:

Output voltage: 5V±5%
Max output current: At least 1A
Output plug polarity: <+> inner
External diameter: 5.5 mm 0.1mm
Internal diameter: 2.1mm 0.1mm


Note:

- When using AC adapter, the power of battery won't be consumed.
- When suddenly stop during measurement (like the plug off from the outlet by carelessness), it must be reinserted the plug into the unit, and restart the measurement.

How to set

1. Start to set, Unit setting:

Press the SET button when power off, 0 or 0.0 will be displayed, then the setting begin. Continue to above step, the unit will be changed when press the MEM button each time. Press the SET button to confirm the unit, then it will enter into the User setting mode.



10

Setting mode

2. Speaking setting

Continue to above step, the screen will display SPE OFF or SPE ON, press the MEM button and it will change between SPE OFF and SPE ON, press the SET button to confirm the option. Following this, the device will enter into the User setting mode.

Remark:
SPE OFF: close the speaking function.
SPE ON: open the speaking function.



3. User setting:

Continue to above step, the screen will display user 1 or user 2, press button MEM, it will be changed between user 1 and user 2, press button SET when you confirm the user, then it will enter into the year setting mode.



4. Year setting:

Continue to above step, the screen will display and flash 20XX, the last digit of the year will increase 1 when press button MEM each time, you could choose from 2020 to 2099. Press button SET when you confirm the year, then it will enter into the month and date setting mode.




11

Setting mode


5. Month and date setting

Continue to above step, the screen will display xxMxxD and xxxx, and keep flashing on month, the digit will increase 1 when press button MEM each time, you could choose from 1 to 12. Press button SET when you confirm the month, then it will set the date. Same as the month setting, each time you press button MEM, the digit will keep changing from 01 to 31. Press button SET when you confirm the date, then it will enter into the time setting mode.



6. Time setting:

Continue to above step, the screen will display xxMxxD and xx:xx, and keep flashing on the digits of hour, the digit will increase 1 when press button MEM each time, you could choose from 0 to 23. Press button SET when you confirm the hour, then the digits of minute start to flash, same as the hour setting, each time you press button MEM the digits will keep changing from 00 to 59. Press button SET when you confirm the minute, then the total setting mode is completed.



Measurement

Pre-measurement

12

Proper use of the unit

- Relax for about five to ten minutes prior to the measurement. Avoid eating, drinking alcohol, smoking, exercising and bathing for 30 minutes before taking a measurement. All these factors will influence the measurement result.

- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Take measurement regularly at the same time of every day, as blood pressure changes even during the day.

Common factors of wrong measurement

- All efforts by the patient to support their arm can increase blood pressure.

- Make sure you are in a comfortable, relax position and do not activate any of the muscles in the measurement arm during measurement. Use a cushion for support if necessary.

- If the arm artery lies lower or higher than the heart, a false reading will be obtained.

Note:

- Only use clinically approved cuffs!

- A loose cuff or a exposed bladder causes false reading.

- With repeated measurements, blood accumulates in the arm which can lead to false reading.

- Consecutive blood pressure measurements should be repeated after 1 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away.

13


Proper Use of the Unit

Fitting the Cuff

- Place the cuff flat on the table with the velcro side down. Pass the end of the cuff through the metal loop so that a circle is formed. The velcro closer will now be facing outwards (ignore this step if the cuff has already been prepared).




- Push the cuff over the left upper arm so that the tube points in the direction of the lower arm.



- Wrap the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow and the rubber tube leaves the cuff on the inner side of the arm.



- Tighten the free end of the cuff and close the cuff by affixing the velcro.



- The cuff should be snug on your upper arm so that you can fit 2 fingers between the cuff and your upper arm. Any piece of clothing restricts the arm and should be removed.



- Secure the cuff with the velcro closer in such a way that it lies comfortably and not too tight. Lay your arm on a table (palm upwards) so that the cuff is at the same height as the heart. Do not bend the tube.



Note:

If it is not possible to fit the cuff to your left arm, it can also be placed on the right. However, all measurements should be made using the same arm.




14

Proper use of the unit

Measuring Procedure:

After the cuff has been appropriately positioned, the measurement can begin as follows:

- Press the START/STOP button, all symbols appear on the display, then the pump begins to inflate the cuff, the rising pressure in the cuff is shown on the display.
- After the suitable pressure has been reached, the pump stops and the pressure gradually falls. The cuff pressure is displayed. In case that the inflation is not sufficient, the device automatically re-inflates to a higher pressure.
- When the device detects the signal, the heart symbol on the display starts to flash.
- When the measurement has been completed, the systolic, diastolic and pulse rate will appear on the display.



Remark: Once finish the measurement, the device will upload the data to the background server automatically with showing two rotating signal bar in the top right corner of LCD (as Pic 1).

After uploading the data successfully, the LCD will show End (as Pic 2) and then the device switches off immediately. If fail to upload the data, LCD will display Err.

Note: The symbol will be displayed along with the reading if the irregular heartbeat is detected during the measurement.

Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (eg. the patient feels unwell) the START/STOP button can be pressed at any time. The device immediately decrease the cuff pressure automatically.

Memory-recall of measurements

This blood pressure monitor automatically stores 2x90 sets measurements value, the oldest record will be replaced by the latest measurement value when more than 90 sets each user.


Read memory record

Press the button MEM when power off, the latest 3 times average value will be shown, press the button MEM again, the last measurement value will be shown, as well as subsequent measurements can be display one after the other by pressing the button MEM each time.



15

About blood pressure



Memory -clear of measurements

If you are sure that you want to permanently remove all stored memories. Press the button SET for 9 times until CL appears when power off, press the START/STOP button, CL will flash for 3 times to clear all the memories. After this press button MEM, [M] and "no" will be shown on the display which mean that no memory in store.

Check IMEI details:

Hold and press the MEM button for 3 seconds when power off, released the button and wait for 5 seconds, then the LCD will display the device IMEI number. And press the START/STOP button to exit.



16

Exceptional Situation

Error Indicators

The following symbol will appear on the display when measuring abnormal.

Symbol	Cause	Correction
E-1	Weak signal or pressure change suddenly	Wrap the cuff properly. Remeasure with correct way.
E-2	External strong disturbance	When near cell phone or other high radiant device, the measurement will be failed. Keep quite and no chatting when measure.
E-3	It appears error during the process of inflating	Wrap the cuff properly. Make sure the air plug is properly inserted into the unit. Remeasure.
E-5	Abnormal blood pressure	Repeat the measurement after relax for 30 mins, if get unusual readings for 3 times, please contact your doctor.
[Battery icon]	Low battery	Connect the USB cable and charger to charge the battery.

Trouble removal

Problem	Check	Cause and solutions
No power	Check the battery power	Charge the battery
No inflation	Whether the plug insert Whether the plug broken or leak	Insert into the air socket tightly Change a new cuff
Err and stop working	Whether move the arm when inflate Check if chatting when measured	Keep the body peaceful Keep quite when measure
Cuff leak	Whether the cuff wrap too loose Whether the cuff broken	Wrap the cuff tightly Change a new cuff

⚠ Please contact the distributor if you can't solve the problem. Do not disassemble the unit by yourself!

19

Warranty Information

Statement

- The intended use: the unit is intended to be used by adults at home or medical center to measure blood pressure and pulse rate from the upper arm.
- The device meets the requirements of IEC 80601-2-30 Part 2 for non-invasive blood pressure monitors.
- Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultatory method, within the limits prescribed by the American National Standard, manual, electronic, or automated sphygmomanometers.
- The risk of patient and user can be lowered to acceptable level.

Warranty Information

- The unit is guaranteed to be free of defects in workmanship and materials under normal use for a period of Five Years from the date listed on the purchase record.
- For repair under this warranty, our authorized service agent must be advised of the fault with the period of the warranty. This warranty covers parts and labor only under normal operations. Any defect resulting from natural causes (e.g., flood, hurricane, etc.) is not within this guarantee. This guarantee does not cover damage incurred by use of the unit not in accordance with the instructions, accidental damage, or being tampered with or serviced by unauthorized service agents.
- Monitor subjected to misuse, abuse, and neglect of this manual's content, non-instructional purposes, unauthorized repair or modifications will be excluded from this warranty.

⚠ The device requires no calibration.

⚠ The device is not repairable and contains no user serviceable parts.

22

EMC Declaration

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity		
Immunity Test	EC 60601-1-2 Test Level	Compliance Level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV 100 kHz repetition frequency	Power supply lines: ±2 kV 100 kHz repetition frequency
Surge IEC 61000-4-5	line(s) to line(s): +0.5kV +1 kV.	line(s) to line(s): +0.5kV +1 kV.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle 70% 25/30 cycles Single phase: at 0 0% 250 cycle (50Hz)	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 250 cycle (50Hz)
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150KHz to 80MHZ: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150KHz to 80MHZ: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz
Radiated RF IEC61000-4-3	10 V/m 80 MHz-2.7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz-2.7 GHz 80 % AM at 1 kHz

NOTE UT is the a.c. mains voltage prior to application of the test level.

25

tellihealth

Arm Blood Pressure Monitor

Shenzhen Union Technology Co., Ltd.
Floor 4-6th of Building D, Jiale Science & Technology Industrial Zone, No. 3, ChuangWei Road, Heshuikou Community, MaTian Street, GuangMing New District, 518106 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

Manufactured for: Telli Health LLC
66 West Flagler Street - Suite 900,
Miami, FL, 33130 USA



Made in China

REV.02

About blood pressure

Signal Strength Indicator:

H means signal is strong, L means signal is weak, _ means there is no signal. These marks will appear after finishing measurement and start to unload data




About blood pressure

Blood pressure is the pressure exerted the arteries.

The systolic blood pressure value represents the blood pressured produced by contraction of the heart muscle.

The diastolic blood pressure value represents the blood pressured produced by relaxation of the heart muscle.



17

Care and Maintenance

Care for the main unit and blood pressure monitor cuff

- Keep the unit in the storage case when not in use.
- Clean the unit with soft dry cloth. Do not use any abrasive or volatile cleaners.
- Never immerse the unit or any component in water.
- Make sure the monitor is off prior to cleaning, a mixture of distilled water and 10 percent bleach could be used.
- Using a spray bottle, moisten a soft cloth towel with the bleach or detergent mix until it is fully saturated. Squeeze any excess moisture from the cloth to avoid any dripping or potential oversaturation of the cuff.
- Wipe all surfaces of the blood pressure monitor cuff thoroughly, making sure to clean the inside and outside of the cuff. Be cautious not to get any moisture in the main unit.
- Using a dry cloth, gently wipe away any excess moisture that may remain on the blood pressure cuff. Lay the cuff flat in an unrolled position and allow the cuff to air dry.

Maintenance

- Do not clean the body and cuff with naphtha, thinner or gasoline, etc.
 - Do not wet the cuff or attempt to clean the cuff with water.
 - Store the unit in a clean and dry location. Do not subject the unit to extreme hot or cold temperature, humidity and direct sunlight.
 - Do not use the device under the influence of electromagnetic interference (nearby cellphones, microwaves, etc.).
- ⚠ We won't be responsible for any quality problem if you don't care for and maintain the product as instructed.

20

EMC Declaration

IEC 60601-1-2: 2014 ME EQUIPMENT and ME SYSTEMS identification, marking and documents for Class B product. Instructions for Use:

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments and so on.

Warning: Don't use near active HF surgical equipment and the RF shielded room of an ME 329mm.n for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the blood pressure monitor, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

If any: A list of all cables and maximum lengths of cables (if applicable), transducers and other ACCESSORIES that are replaceable by the RESPONSIBLE ORGANIZATION and that are likely to affect compliance of the ME EQUIPMENT or ME SYSTEM with the requirements of Clause 7 (EMISSIONS) and Clause 8 (IMMUNITY). ACCESSORIES may be specified either generically (e.g., shielded cable, load impedance) or specifically (e.g., by MANUFACTURER and EQUIPMENT OR TYPE REFERENCE).

If any: The performance of the ME EQUIPMENT or ME SYSTEM that was determined to be ESSENTIAL PERFORMANCE and a description of what the OPERATOR can expect if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES (the defined term "ESSENTIAL PERFORMANCE" need not be used).

23

EMC Declaration


Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity																					
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE REPORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	DISTANCE (m) TEST LEVEL (V/m)															
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27															
							450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28								
														710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9	
																					745
							810	870	800-960	Pulse modulation 18 Hz	2	0.3	28								
														1720	1845	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
							1970	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3								
														5240	5500	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
							5785														


26

About blood pressure

- According to the blood pressure classification by the WHO/ISH,
- SYS lower than 100mmHg (13.3kPa) is considered as hypotension.



Blood pressure type



18

Specification

Description	Automatic upper arm blood pressure monitor	
Display	LCD digital display	
Measuring principle	Oscillometric method	
Measuring localization	Upper arm	
Measurement range	Pressure	0~299mmHg
	Pulse	40-199 pulses/min
Accuracy	Pressure	±3mmHg
	Pulse	±5% of reading
LCD indication	Pressure	3 digits display of mmHg
	Pulse	3 digits display
Memory function	Symbol	Memory/Heartbeat/Low battery
	Memory function	2x90 sets memory of measurement values
Power source	1pc 3.7V 1000mAh rechargeable lithium battery	
Automatic power off	in 3 minutes	
Main unit weight	Approx. 219g (batteries not included)	
Main unit size	132mmx100mmx45mm	
Main unit lifetime	10,000 times under normal use	
Battery life	Could be used for about 70 times when charge is full	
Accessories	Cuff, instruction manual, lithium battery, USB cable, Wall charger	
Operating environment	Temperature	5°C-40°C
	Humidity	15%~93%RH
	Air pressure	86kPa - 106kPa
Storage environment	Air pressure: 86kPa ~ 106kPa; Temperature: -20°C ~ 55°C; Humidity: 10% ~ 93%RH; avoid crash, direct sun or rain during transportation	
Expected service life	Five years	
Software version	UA1.0	

21

EMC Declaration

Technical Description

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.

2. Guidance and manufacturer's declaration – electromagnetic emissions and Immunity.

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions	
Emissions Test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Compliance

FCC Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC RF Exposure Information and Statement

When carrying the product or using it while worn on your body, either use an approved accessory such as a holster or otherwise maintain a distance of 5 mm from the body to ensure compliance with RF exposure requirements. Use of other accessories may not ensure compliance with FCC RF exposure guidelines.

27